

Accessories

- ▶ *Sample holders - standard models for typical biological samples or custom designed for unique applications.*
- ▶ *Stage mounting - standard mounting plates for microscopes as well as custom mounting for vertical and horizontal orientations.*
- ▶ *Work directly with Mad City Labs' design engineers to create the accessories needed for your experimental setup at minimal additional cost.*
- ▶ *Avoid in-house machine shop delays*

Microscope Mounting Adapters

Mount nanopositioners directly onto microscopes without a coarse positioning stage



Olympus IX inverted microscope series mount for the Nano-LP Series and Nano-Bio Series




Olympus IX inverted microscope series mount for the Nano-H Series



Olympus IX series, Zeiss Axiovert, and Leica DMI inverted microscope mount for the Nano-LP Series, Nano-Bio Series, and 2-axis Nano-T Series.

Biological Sample Holders

**Compatible with Nano-BioS
and Nano-LPS Series stages**



Re-entrant 35mm petri dish holder

Re-entrant 3 inch (75mm) slide holder

Compatible with 2.6" x 2.6" aperture stages




Re-entrant petri dish holder



Re-entrant coverslip holder



Re-entrant 2 inch slide or coverslip holder



Top mount 3 inch slide holder

Compatible with the Nano-Z500

Note: These designs can be modified to fit other stages



3 inch slide holder



Multimell plate holder



Coverslip adapter



35mm petri dish adapter