

Accessories

- ▶ *Sample holders - standard models for typical biological samples or custom designed for unique applications.*
- ▶ *Stage mounting - standard mounting plates for microscopes as well as custom mounting for vertical and horizontal orientations.*
- ▶ *Work directly with Mad City Labs' design engineers to create the accessories needed for your experimental setup at minimal additional cost.*
- ▶ *Avoid in-house machine shop delays*

Microscope Mounting Adapters

Mount nanopositioners directly onto microscopes without a coarse positioning stage



Olympus IX inverted microscope series mount for the Nano-LP Series and Nano-Bio Series




Olympus IX inverted microscope series mount for the Nano-H Series



Olympus IX series, Zeiss Axiovert, and Leica DMI inverted microscope mount for the Nano-LP Series, Nano-Bio Series, and 2-axis Nano-T Series.

Biological Sample Holders

Compatible with 2.6" x 2.6" aperture stages:



Compatible with Nano-BioS and Nano-LPS Series stages:



Compatible with the Nano-Z500:

Note: These designs can be modified to fit other stages

