



UHV Nanopositioning

December 2015

Capacitive sensors perform poorly in vacuum, with performance degrading up to 60% per meter of cable length. Mad City Labs **PicoQ[®] sensor performance does not degrade in vacuum**. Our extensive and practical knowledgebase, combined with proprietary **PicoQ[®] sensor technology**, uniquely places us to provide a standard product line of **UHV nanopositioners and micropositioners** and develop customized solutions with the highest performance characteristics.

- **lowest UHV prices**
- best UHV design
- highest resolution available in UHV
- **PicoQ[®] sensor technology**
- Unlike capacitive sensors, the **performance of PicoQ sensors are not degraded in UHV**
- UHV engineers with over 50 years combined UHV experience
- Standard nanopositioning systems can be customized



3 axis custom UHV system.


Applications

- X-ray, VUV, and optical microscopy and spectroscopy
- Surface metrology
- UHV atomic scale microscopy
- Interferometry
- Particle physics




Custom three axis tip/tilt nanopositioning solution for UHV 10 meter atomic interferometer.

Mad City Labs has a standard line of **UHV compatible nanopositioners**. These nanopositioners have a central aperture for beam lines and optical transmission with 50 to 200 μ m range of motion on two or three axes.



Nano-UHV200 system with 200 μ m range of motion in XYZ, closed loop control, and large central aperture. Bakeable to 100 C.

Product Catalog



Website



Unsubscribe

You are receiving this message because you have expressed interest in Mad City Labs products. If you do not wish to receive future newsletters, please click on the red unsubscribe link above or reply to this message with "Unsubscribe from Mad City Labs Newsletters" in the subject.

2524 Todd Drive
Madison, WI 53713 USA

USA +1 (608) 298-0855
Europe +41 (0)58 269 8017

information@madcitylabs.com
www.madcitylabs.com

© 2015